

BaseScope™ Probes for MAGEA Gene Family Detection

Discern highly homologous sequences in intact tissue with morphological context

Featured publications using the BaseScope™ Assay for short or highly homologous RNA targets:

Functional ectopic neuritogenesis by retinal rod bipolar cells is regulated by miR-125b-5p during retinal remodeling in RCS rats.
Fu et. al. Sci Rep. 2017; 7(1): 1011

Rare progerin-expressing preadipocytes and adipocytes contribute to tissue depletion over time.
Revechon et. al. Sci Rep. 2017; 7(1): 4405

Involvement of DHH and GLI1 in adrenocortical autograft regeneration in rats.
Takizawa et. al. Sci Rep. 2018; 8(1): 14542

Lariat intronic RNAs in the cytoplasm of vertebrate cells.
Talhouarne and Gall. Proc Natl Acad Sci U S A. 2018; 115(34): E7970

See more publications at acdbio.com/publications

MAGEA Gene Family Detection

Ongoing clinical trials of cancer vaccines and adoptive cell therapies target members of the melanoma-associated antigen (MAGE) family, which are highly prevalent in tumors. However, their highly homologous nature, with >90% sequence identity between the *MAGE-A* genes, poses significant challenges to their specific detection in tumors. Due to antibody cross-reactivity there is limited capability of protein-specific assays to detect and distinguish between the various MAGE-A antigens. Here, the BaseScope™ RNA *in situ* hybridization (ISH) assay was developed to discern *MAGEA1*, *MAGEA3*, *MAGEA4*, and *MAGEA10* expression in cell lines and tumor tissues.

The BaseScope™ Assay is similar to the RNAscope® ISH technology. Both achieve single-molecule RNA detection using paired oligo ("ZZ") probes to amplify signal without non-specific background. However, the BaseScope™ probe design and advanced signal amplification system enable single-molecule RNA detection with as few as 1 ZZ probe and it can differentially discern between highly homologous sequences, such as the *MAGE-A* gene family members, in intact fixed tissue.




FIGURE 1. Detection of individual *MAGE-A* genes in human lung cancer tissue.




FIGURE 2. Specific detection of individual *MAGE-A* genes in human cancer cell lines.

Gene	H460* (Lung Cancer)	MBA-MB-468 (Breast Cancer)	A431 (Skin Cancer)	SK-MEL-2 (Skin Cancer)	SK-OV-3* (Ovarian Cancer)
<i>MAGEA1</i>	+++++	-	-	-	-
<i>MAGEA3</i>	+	+++++	-	++	-
<i>MAGEA4</i>	-	-	+++++	-	-
<i>MAGEA10</i>	-	+	-	+++++	-

TABLE 1. Summary of *MAGEA* expression pattern in cells lines as detected by the BaseScope™ Assay.

*H460 and SK-OV-3 cell line expression patterns of *MAGEA1*, *MAGEA3*, *MAGEA4*, and *MAGEA10* confirmed by GeneCards







FIGURE 3. Detection of individual MAGE-A genes in human melanoma.







FIGURE 4. Detection of individual MAGE-A genes in human head & neck cancer.




FIGURE 5. BaseScope™ Assay workflow.

Probes	Catalog #	Comments
BA-Hs-MAGEA1-2zz-st	715201	Detects <i>MAGEA1</i>
BA-Hs-MAGEA3-2zz-st	715211	Detects <i>MAGEA3</i>
BA-Hs-MAGEA4-2zz-st	715221	Detects <i>MAGEA4</i>
BA-Hs-MAGEA10-2zz-st	715231	Detects <i>MAGEA10</i>

Assay	Catalog #
BaseScope™ LS Reagent Kit	323600

TABLE 2. BaseScope™ probes and reagents used for *MAGE-A* detection.

Need to know what MAGEA family members your tissues are expressing? Outsource your study to the experts!

- Utilize ACD's expertise in the BaseScope MAGEA assay for recommendations on a controlled study design
- Rely on ACD's scientific experts for accurate data interpretation
- Results are delivered in 4 weeks on average
- Have confidence in our quality procedures and optimization protocols

"This data puts us in a perfect position to act immediately so thank you for all of your effort on this and its successful execution."

– Scott Rowlinson, Aeglea Bio



For more information, visit us online at www.acdbio.com

