

Setting new standards in biological safety cabinets





Thermo Scientific biological safety cabinets are designed to maximize sample protection and user safety, with an emphasis on containment, comfort and convenience, all day, every day. For more than 5 decades we have been a pioneer and innovator of clean air technologies. Our global vision and experience have enabled us to continually provide the best application based solutions ideally suited for your requirements and working environment, with compliance to all major worldwide safety standards.



Precise self adjusting protection

Autonomous safety features maintain proper air flow and monitor critical conditions in real time without user interaction. This supports valuable research enabling our customers to focus where needed.



Simplified operation and ergonomics

Readily understandable graphical display provides the information you need, while electronic sash with fully hinged opening enables easy access for supplies and equipment, and simplifies easy cleaning and disinfection



Enhanced information management

Electronically collect key operational and performance data on the secure Thermo Scientific Cloud and remotely monitor from anywhere, to evaluate your processes and accelerate your discoveries


Simplified operation and ergonomics



Readily understandable graphical display provides the information you need, while electronic sash with fully hinged opening enables easy access for supplies and equipment, and enables easy cleaning and disinfection.

Helping ensure that the BSC is managing safety conditions autonomously, and is simple and intuitive to use allows our customers to fully focus on demanding applications and precision science.

Intuitive | Simple | Convenient



Fully opening front-hinged window design enables easy cleaning and disinfection with full access to the interior cabinet space.

Full opening also allows convenient posture for thorough cleaning and visible reach to all cabinet surfaces including rear wall.

Most resistant to scratches and chemical attack in the user zone. All surfaces inside the HeraSafe 2030i are from 316L stainless steel The finish is "brushed and sanded" which both reduces glare and offers better scratch resistance.

Smart technology delivering optimal containment and energy efficiency

Autonomous safety features insure continuous accuracy that never takes a day off – providing optimal protection from contamination by compensating to match the way you work and lab environment

Sash Opening 200/250mm | Smart Flow+ and Dave+ | Password Protection | UV Protection Energy | Efficient Design |




Electronic front sash window

with push button simplicity the front window automatically positions itself at the position mandated by local regulation to assure proper working conditions. Sloped front window allows fatigue free working position and is aerosol tight sealed when closed

200mm sash opening (Standard) provides optimum balance of comfort plus efficiency.

250mm sash opening (optional) provides additional convenience for quick and safe consumables movement.





SmartFlow Plus

patented proprietary design automatically compensates unique energy saving dual DC blower motors to insure the proper airflow rates no matter the loading condition of the protective HEPA filters. Unlike typical BSC's which rarely perform at the same levels day after day due to HEPA loading, and potential disruptions to airflow such as paper covering vent openings.

How does this work?

Independent intelligent DC motors for downflow and inflow monitor and control fan speed in real time to maintain balanced airflow during filter loading or temporary obstruction.

Exclusive to our BSCs, independent supply and exhaust blowers automatically balance both the downflow and inflow velocities to maintain both operator and sample protection.

Dave Plus


our unique digital airflow verification system validates protection by monitoring airflow velocities with independent pressure sensors which detect and alert to any changes across the exhaust and downflow plenums.

How does this work?

By sensing pressure in the throat of the blowers, DAVePlus monitors airflow directly for the inflow and downflow in the Herasafe 2030i to assure both personal and sample protection.


Controlled Passcode Access

for convenience the BSC can be configured via the user interface to operate in "open mode" or "secure mode" depending upon working preferences and lab environment. In secure mode a selectable passcode must be entered to enable operation of the BSC. In this way different units utilized by different personnel for specific applications, may be controlled by providing a unique passcode to each, or the same code to all.



Self adjusting germicidal UV light intensity

for those who routinely employ UV germical light to their work surface at the end of the work day, we provide a time adjustable functionality to auto shut off the light at the desired interval to preserve bulb life. Our exclusive self adjusting exposure capability insures reproducible germicidal effectiveness with every use, by assuring performance at consistent levels to protect your important work.




Energy efficient design

- Energy conserving dual DC motor design, generates less heat than older AC designs with double the usage lifetime and digitally controllable for precision performance
- Stand-by mode automatically reduces blower speed to protect samples with reduced energy usage until you are ready to resume work

Optimizing total cost of operation

Thermo Scientific SmartFlow airflow compensation and balanced filters have led to safer BSCs that are easier to use and significantly less costly to operate.**

The more you use it, the more you save!



Costs, published in American Biotechnology Laboratory, September 2008.

^{**} Compared to previous models.

^{***} Comparison based on 15 year life of 1.2 m wide biological safety cabinets and 2011 average costs of energy and cooling in the EuroZone.

"Traditional BSC" estimates based on published product specifications for readily available BSCs using AC motor technology.

Innovative productivity enhancements simplify use, comfort and cleaning



Full color touch screen user interface – safety you can see

Our advanced yet simple to operate, best in class graphic user interface is customer tested to provide easy navigation and convenient illustration of airflow safety and usage data to address any lingering user concerns about the cabinet functioning properly to product them and their valuable samples. Performance details such as user interaction logs, errors and event logs are easily identified to drive continuous process improvement and reproducibility of results.



Standard smart ports

provided as standard with every cabinet, Smart Ports provide clean and safe approach to incorporating vacuum tubing and cables through the side wall of the BSC, for improved organization and work efficiency. The side wall plenum is under negative pressure preventing possibility of leakage in or out of the BSC during operation.

Customisation

options are available for the Herasafe 2030i. We have vast experience in designing custom adaptations to our BSC's. We can design from scratch or provide frequently requested solutions.



Large interior work space/premium 316 non-corrosive stainless steel

conveniently large interior workspace provides ample room for incorporation of instruments and process optimization devices including robotic arms. We utilize a "brushed and sanded" finish to provide maximum scratch resistance and non-reflective properties.

Padded armrests

adjustable ergonomic armrests are provided with each unit to assure proper arm support to reduce strain during prolonged work periods while extending reach within the cabinet. Removable gel filled comfort pads are available as an option to enhance the user experience.

Interior LED lights

visibility within the cabinet is enhanced with the use of bright, energy efficient lighting, eliminating the potential for fatigue often experienced with use of standard fluorescent lighting

Exclusive one touch start /stop function

exclusive one touch auto start and stop function activates or deactivate all cabinet components efficiently, and ready to use without lengthy wait times.

Smart Clean hinged front window

our exclusive fully opening front hinged window design enables easy cleaning and full access to the cabinet interior work space for thorough cleaning/disinfection and incorporation of ancillary work tools and equipment as may be desired.



Enhanced information management with the convenience of connectivity



Most of the enhanced information provided by the touch screen interface, can be remotely monitored and captured by web based or mobile app for subsequent analysis, eliminating the need for paper notebooks

Standard USB data outputs – enable export of data to external data systems if desired. A powered connection is included in the front panel to enable convenient accessory charging and accommodate future software update uploads

Convenient event and alert messaging logs -

accessible in the user interface provide useful and important usage data and alert information for process optimization and maintenance

Embedded connectivity – ready for plug and play into the Thermo Fisher secure cloud as desired. Enables simplified information sharing with colleagues and enhances productivity through performance monitoring from anywhere, anytime improving process control, reporting and asset management, with use of the exclusive free Thermo Scientific mobile app

Future proof – once you have incorporated the Thermo Scientific Cloud you are positioned to benefit from continuous improvement of data management, as your capabilities with expand as the Thermo Scientific cloud evolves to include new functionality over time (active Q4 2018)

Herasafe 2030i Ordering Information			EN 12469				NSF/ANSI 4	9	
Size/Internal width		m (ft)	0.9 (3)	1.2 (4)	1.5 (5)	1.8 (6)	1.2 (4)	1.8 (6)	
200mm Sash with Cross Beam UV-C in side v	valls		51032328	51032329	51032330	51032331			
200mm Sash – no options included** (250mm/10" sash can be optioned)			51032324	51032325	51032326	51032327			
250mm (10") Sash with Cross Beam UV-C in s single piece work tray and prep for 3x rear wal							51032334	51032335	
250mm (10") Sash, indented single piece work	tray						51032332	51032333	
Window operating aperture 10" (250mm)		Option	51901247**				Sta	ndard	
Dimensions									
External Dimension (Cabinet only) d	epth	mm (in.)	800 (31.5)	800 (31.5)	800 (31.5)	800 (31.5)	800 (31.5)	800 (31.5)	
	vidth	mm (in.)	1000 (39.4)	1300 (51.2)	1600 (63.0)	1900 (74.8)	1300 (51.2)	1900 (74.8)	
h	eight	mm (in.)	1536 (60.5)	1536 (60.5)	1536 (60.5)	1536 (60.5)	1536 (60.5)	1536 (60.5)	
	lepth	mm (in.)	465 (18.3)	465 (18.3)	465 (18.3)	465 (18.3)	465 (18.3)	465 (18.3)	
(net, without air intake grill)	vidth	mm (in.)	900 (35.4)	1200 (47.2)	1500 (59.1)	1800 (70.9)	1200 (47.2)	1800 (70.9)	
Shipping Dimension d	lepth	mm (in.)	925 (36.4)	925 (36.4)	925 (36.4)	925 (36.4)	925 (36.4)	925 (36.4)	
	 vidth	mm (in.)	1110 (43.7)	1410 (55.5)	1710 (67.3)	2010 (79.1)	1410 (55.5)	2010 (79.1)	
h	eight	mm (in.)	1710 (67.3)	1710 (67.3)	1710 (67.3)	1710 (67.3)	1710 (67.3)	1710 (67.3)	
Working height of front window		mm (in.)	, ,	, ,	00 (8)		, ,	60 (10)	
Maximum window opening (for cleaning, loading/unloading)		mm (in.)	773 (30.4)	773 (30.4)	773 (30.4)	773 (30.4)	773 (30.4)	773 (30.4)	
Work surface height		mm (in.)		750-950	(29.5-37.4)		750-950	(29.5-37.4)	
Weight									
Net weight		kg (lbs)	170 (375)	200 (441)	230 (507)	280 (617)	200 (441)	280 (617)	
Shipping weight		kg (lbs)	190 (419)	225 (496)	260 (573)	315 694)	225 (496)	315 694)	
Max load capacity of one-piece work surface		kg (lbs)	50 (110)	50 (110)	50 (110)	50 (110)	50 (110)	50 (110)	
Max load capacity of single segmented work s (per 30cm segment)	surfaces	kg (lbs)	25 (55)	25 (55)	25 (55)	25 (55)	25 (55)	25 (55)	
Electrical Data									
Voltage		V	230 120						
Hertz		Hz			60				
Protection class			1/IP20						
Filter Specification									
Specification				H1	4 HEPA EN18	22, 99.995%	MPPS		
Certifications									
Standard							ANSI 49		
Agency certified			TÜV *	TÜV	TÜV *	TÜV	NSF		
Certification LNE (France only)		Option		LNE			LNE		
Ergonomics and Utilities									
Lighting intensity		lux (fc)	1000 (93)	929 (86)	1500 (139)	1418 (132)	850 (79)	1300 (121)	
Interior outlets		pcs.	2 rea	ar wall		ar wall	2 rear wall	4 rear wall	
Service valve ports		pcs.	2 on each side + Smart Port				2 on each si	de + Smart Po	
Slope of front window		degree °			10			10	
Herasafe 2030i Specifications									
Noise level (measured according to EN 12469) 200mm opening		dB (A)		57		57			
Noise level (EN units measured according to ENSF units measured according to NSF/ANSI 4 250mm opening				57		58	62	62	
Power consumption (in normal operation, fans lights on) 200mm opening	on,	W/h	133	153	217	251			
Power consumption (in normal operation, Fans on, lights on) 250mm opening			155	175	240	300	180	316	
Power consumption (in stand-by mode)		W/h	55	55	80	80	55	80	
Normal heat output in operation mode (non-vented) 200mm opening		W/h	133	153	217	251			
Normal heat output in operation mode (non-vented) 250mm opening			155	175	240	300	180	316	
Stand-by heat output (lights off, non-vented)		W/h	55	55	80	80	55	80	
* TÜV/CC quailable for madele with 200mm or				'					

^{*} TÜV GS available for models with 200mm access opening

Herasafe 2030i Ordering Information		EN 12469				NSF/ANSI 49	
Size/Internal width	m (ft)	0.9 (3)	1.2 (4)	1.5 (5)	1.8 (6)	1.2 (4)	1.8 (6)
Exhaust volume of cabinet 200mm opening	m3/h (cf/m)	296 (174)	395 (232)	493 (290)	592 (348)		
Exhaust volume of cabinet 250mm opening	m3/h (cf/m)	370 (218)	494 (291)	617 (363)	741 (436)	582 (342)	872 (513)
Exhaust volume, thimble ducting 200mm opening	m3/h (cf/m)	385 (226)	513 (302)	641 (377)	770 (453)		
Exhaust volume, thimble ducting 250mm opening	m3/h (cf/m)	481 (283)	642 (378)	802 (472)	963 (567)	756 (445)	1134 (667)
Recentacle Ontions							

Country configurations are available to ensure correct line cord and internal receptacles are shipped with the cabinet. Standard version for EN12469 cabinets is Schuko, and for NSF cabinets is type B

0	0 - 1'	E4000000	
Country version Switzerland	Option	51900300	
Country version UK	Option	51900303	
Country version Italy	Option	51900306	
Country version FR/B/CZ/SK/PL	Option	51900921	
Country version Australia	Option	51900449	
Country version Denmark	Option	51900481	
Country version China	Option	51900900	
Country version Israel	Option	51901255	
Country version India	Option	51901261	
Country version Argentina			51900308 (230V)
Country version Chile			51901259 (230V)
Country version Brazil			51901260 (230V)
0. 10.0			

Stand Options

We recommend ordering a base stand to ensure a robust and stable installation of the BSC. All of our base stands have been tested with our units for stability and functionality. Choose between a fixed height, manually adjustable working height, or electrically adjustable working height. Please consider the total height of the room when choosing the right working height. Footrest and castor options are available for added convenience.

Fixed height stand charcoal grey Height of 680 - 880 mm	Acc.	50155665	50155666	50155667	50155668	50155666	50155668
Adjustable height stand charcoal grey Height of 680 - 880 mm	Acc.	50155688	50155689	50155690	50155691	50155689	50155691
Electric Adjustable Stand grey-white (RAL 9002)	Acc.	50117449	50109314	50117450	50109325	50109314	50109325
Footrest bar	Acc.	50051982	50051983	50051984	50051985	50051983	50051985
Castor-Set for Stand grey-white (4 pcs.) (RAL9002) Provides work height of 850 mm	Acc.		5010)9977		501	09977

Work Surface Options

You have the option to fully customize your new BSC to accommodate any work task. The standard EN cabinet includes our stainless steel segmented worktray for easy removal, autoclaving, and cleaning under the work surface. The standard NSF cabinet includes our stainless steel single piece indented worktray. To customize your work area, order the part number to remove the standard worktray and configure the cabinet as desired. All internal surfaces and worktrays are 316 SS.

Without table tops (negative price)		Option	51900651	51900652	51900653	51900654	51900652	51900654	
Flat worksurface module (150 mm)	KAM15	Acc.		500	73666		500	73666	
Flat worksurface module (300 mm)	KAM30	Acc.		500	500	72685			
Flat worksurface module (600 mm)	KAM60	Acc.		500	500	50073667			
Flat worksurface module (900 mm)	KAM90	Acc.		500	50076398				
Lowered worksurface module (600 mm)	KAML60	Acc.	50073668				500	50073668	
Lowered worksurface module (900 mm)	KAML90	Acc.	50073669				500	50073669	
Lowered worksurface module (1200 mm)	KAML120	Acc.	No 50073670				500	073670	
Lowered worksurface module (1500 mm)	KAML150	Acc.	No 50073671		No	50073671			
Lowered worksurface module (1800 mm)	KAML180	Acc.		No	<u>'</u>	50078750	No	50078750	
Vibration dampened surface module with granite plate (300 mm)	KAMW30	Acc.		50073757				50073757	
Vibration dampened surface module with granite plate (600 mm)	KAMW60	Acc.		50076162				50076162	
3 liter mini-sink table top module	KAMS30	Acc.		500	78875		50078875		
Replacement armrests KAS 30 (Set of 2 armrests 316 SS)	KAS 30		50073663				500	73663	
Padded armrest 316 SS (Set of 2 including gel inserts)		Acc.		50156438				56438	
Replacement gel pads for padded armrest	(2 pcs.)	Acc.		501	56437		501	56437	

Herasate 2030i Ordering Information	mc /f1\	EN 12469	1.0 (4)	1 5 (5)	1.0.(0)	NSF/ANSI 4		
Size/Internal width	m (ft)	0.9 (3)	1.2 (4)	1.5 (5)	1.8 (6)	1.2 (4)	1.8 (6)	
Additional Filter Options								
Additional exhaust HEPA filters (AEF) can be fitted to filters are not a requirement for recirculating BSC's, of carbon filters can be fitted on top of the BSC, but sh	even in the UK, v	where BS5726	6 has been su	perceded by E	EN12469, sinc	e 2000. Addit	ional exhaust	
Additional exhaust HEPA Filter, retro-fit KAEF	Acc.	50078857	50073772 50078856		50078856			
Additional exhaust Carbon Filter, retro-fit KAEFC	Acc.	50078783	50073773	50073773 50078784				
Grid filter set, G1 acc. To EN 779 (5 pcs)		50076846		50076921		50076846	50076921	
JV Options							1	
The UV control function is designed to preserve bulk measures and tracks UV bulb use ensuring that the Additionally by monitoring UV bulb age, the system a	bulbs will not re	quire excessiv	ve replaceme	nt due to prem				
Cross beam UV-C light option in side walls ncludes protective cover								
ncluded in models: 51032328, 51032329, 51032330, 51032331	Standard	Sta	ndard (Cross	Beam - in side	walls)			
Available for models: 51032324, 51032325, 51032326, 51032327	Option	5190	00642 (Cross	Beam - in side	e walls)			
ncluded in models: 51032328, 51032329, 51032330, 51032331	Standard					Standard (cross beam - in side wa		
Available for models: 51032324, 51032325, 51032326, 51032327	Option					51900642 (cross beam - in side wa		
Media Valves								
Our standard BSCs feature a dual-side wall negative also offer factory-installed piping options on our pre								
Rear Wall Media Connections								
Butility connectors prep (rear wall) Select rear wall service taps additionally	Option							
wailable for models: 51032324, 51032325, i1032326, 51032327	Option		519012	251 (G3/8")				
ncluded in models: 51032328, 51032329, i1032330, 51032331	Standard					Included		
Available for models: 51032324, 51032325, 51032326, 51032327	Option					51901252 (NPT3/8")		
Rear Wall Service Taps/Media Valves (supplies a	re connected a	above BSC)						
Combustible gas tap for rear wall	Acc.		500	068235		500	068235	
Ion-combustible gas tap for rear wall	Acc.		500	068234		50068234		
Vater tap for rear wall	Acc.		500	068233		50068233		
acuum tap for rear wall	Acc.		500	068232		50068232		
Side Wall Service Taps/Media Valves/Lead Thro	ughs							
N 1 171 1 6 11 11	Acc.		500	046015		50046015		
Combustible gas tap for side wall			50045959			50045959		
	Acc.		500	045959		500	045959	
lon-combustible gas tap for side wall	Acc.)45959)44679)45959)44679	
Ion-combustible gas tap for side wall Vater tap for side wall			500			500		
Vacuum tap for side wall	Acc.		500 500	044679		500 500	044679	
Vater tap for side wall Vacuum tap for side wall, hose dia. 10-13mm	Acc.		500 500 500	044679 044678		500 500 500)44679)44678	
Vater tap for side wall Vacuum tap for side wall Hose connection for side wall, hose dia. 10-13mm Hose connection for side wall, hose dia. 3mm	Acc. Acc. Acc.		500 500 500	044679 044678 059017		500 500 500	044679 044678 059017	
Vater tap for side wall Vater tap for side wall Vacuum tap for side wall Hose connection for side wall, hose dia. 10-13mm Hose connection for side wall, hose dia. 3mm Hose connection for side wall, hose dia. 6/3mm	Acc. Acc. Acc. Acc. Acc.	nnection	500 500 500	044679 044678 059017 076408		500 500 500	044679 044678 059017 076408	
Non-combustible gas tap for side wall Vater tap for side wall Vacuum tap for side wall Hose connection for side wall, hose dia. 10-13mm Hose connection for side wall, hose dia. 3mm Hose connection for side wall, hose dia. 6/3mm Solenoid Valve (should be quoted with combusti	Acc. Acc. Acc. Acc. Acc.	nnection	500 500 500 500 500	044679 044678 059017 076408		500 500 500 500 500	044679 044678 059017 076408	
Vater tap for side wall Vacuum tap for side wall Hose connection for side wall, hose dia. 10-13mm Hose connection for side wall, hose dia. 3mm Hose connection for side wall, hose dia. 6/3mm Solenoid Valve (should be quoted with combusti	Acc. Acc. Acc. Acc. Acc. blegas G3/8 cc	nnection	500 500 500 500 500	044679 044678 059017 076408 066668		500 500 500 500 500	044679 044678 059017 076408 066668	
Non-combustible gas tap for side wall Vater tap for side wall Vacuum tap for side wall Hose connection for side wall, hose dia. 10-13mm Hose connection for side wall, hose dia. 3mm Hose connection for side wall, hose dia. 6/3mm Solenoid Valve (should be quoted with combusti Solenoid valve for combustible gas tap Exhaust Duct Connection Kits HEPA filters do a great job of trapping particles and nicrobiological studies. In these applications, it is im	Acc. Acc. Acc. Acc. Acc. blegas G3/8 cc Acc. biohazards, but	they do not tr	500 500 500 500 500 7ap gases and duct transition	044679 044678 059017 076408 066668 074905 d toxic chemica n to route exha	ust air out of t	500 500 500 500 500	044679 044678 059017 076408 066668 074905	
Non-combustible gas tap for side wall Nater tap for side wall Accuum tap for side wall Hose connection for side wall, hose dia. 10-13mm Hose connection for side wall, hose dia. 3mm Hose connection for side wall, hose dia. 6/3mm Solenoid Valve (should be quoted with combusting Solenoid valve for combustible gas tape Exhaust Duct Connection Kits HEPA filters do a great job of trapping particles and microbiological studies. In these applications, it is impervironment. We offer two different connections that Thimble Duct Exhaust Transition KTD	Acc. Acc. Acc. Acc. Acc. blegas G3/8 cc Acc. biohazards, but	they do not tr	500 500 500 500 500 rap gases and duct transition erent ways, de	044679 044678 059017 076408 066668 074905 d toxic chemica n to route exha	ust air out of t	500 500 500 500 500	044679 044678 059017 076408 066668 074905	
Combustible gas tap for side wall Non-combustible gas tap for side wall Water tap for side wall Vacuum tap for side wall Hose connection for side wall, hose dia. 10-13mm Hose connection for side wall, hose dia. 3mm Hose connection for side wall, hose dia. 6/3mm Golenoid Valve (should be quoted with combustice of the combustice	Acc. Acc. Acc. Acc. blegas G3/8 cc Acc. biohazards, but portant to select can be installed	they do not tr t an exhaust d in three diffe	500 500 500 500 500 500 crap gases and duct transition erent ways, de 500	044679 044678 059017 076408 066668 074905 d toxic chemica in to route exha	ust air out of t facility requir	500 500 500 500 500 600 e used as an athe laboratory ements.	044679 044678 059017 076408 066668 074905 adjunct to to the extern	

EN 12469

NSF/ANSI 49

Herasafe 2030i Ordering Information

Herasafe 2030i Ordering Information			EN 12469		NSF/ANSI 49			
Size/Internal width		m (ft)	0.9 (3)	1.2 (4)	1.5 (5)	1.8 (6)	1.2 (4)	1.8 (6)
User Convenience Options								
We offer unique customer installed access	ories to help	organize you	r workspace a	and free up va	luable space i	nside and aro	und the BSC.	
Hanging bar with set of 6 hooks	IVB 6 6	Acc.	50073943	50073943 No			No	
Hanging bar with set of 12 hooks	IVB 10 12	Acc.	No	50073944		50073944		
Utility basket to be used with hanging bar	KB	Acc.		50078900				78900
2 utility baskets incl. additional	KBB	Acc.		50078901				78901

Herasafe 2030i Ordering Information		
Additonal topboxes for fumigation	Part number	m (ft)
	KSF9	Acc.
KSF frame fixed to top of each recirc cabinet to accept the relevant topbox for fumigation. Brackets are fitted to side of frame to allow the use of speed clamps. Height = 50mm	KSF12/15	Acc.
nted to didd of marrie to allow the add of oppose stamps. Holgrit – softim	KSF18	Acc.
KSNB topbox for use with Hydrogen Peroxide fumigation.Speed clamps allow for fixture	KSNB9	Acc.
to a KSF frame. The 35mm hole will accept a Steris OR Bioquell outlet camlock fitting.	KSNB-12/15	Acc.
Height = 110mm with gasket.	KSNB-18	Acc.
KS-EVBOX for use with formonia system and formaldehyde. Speed clamps allow it to be clamped to	KS-EVBOX-9	Acc.
KSF frame, while the 2 threaded holes accept the ABEKP3 charcol filters. This box can only be used if	KS-EVBOX-12/15	Acc.
the formaldehyde is being neutralized with ammonia after fumigation. Height = 110mm with gasket.	KS-EVBOX-18	Acc.
KSDB for use with formaldehyde fumigation only. Speed clamps allow for fixture to a KSP frame. This	KS9DB160	Acc.
top box is supplied with a 3m flexi trunking which allows you to vent formaldehyde out into another exhausted application to purge fumes from cabinet.	KS12DB160	Acc.
Side exit only. Height = 190mm with gasket.	KS18DB160	Acc.
KSP(FD) For use with formaldehyde fumigation. Fixes to a KSF frame. Supplied with a	KSP(FD)-9	Acc.
3m flexi trunking which allows venting of fumes to an exhausted application. Top exit	KSP(FD)-12/15	Acc.
only. Height = 50mm (Frame).	KSP(FD)-18	Acc.
ABEKP3 charcoal filters screw into KS-EVBOX. Pack comprises of 3 filters.	ABEKP3	Acc.

Not all BSCs are manufactured to the same safety standards. Our design philosophy is grounded in compliance with the recognized worldwide safety standards, because your work cannot compromise on safety. Our biological safety cabinets are subjected to all the following certifications:



hanging bar

EN 12469

Our Class II biological safety cabinets are independently tested and certified to EN 12469 standard for Class II biological safety cabinets by TÜV Nord, the leading testing agency in Europe.



Our Class II biological safety cabinets are independently tested and certified to EN 12469 by LNE, the leading testing agency in France.



Our Class II biological safety cabinets are independently tested and certified to NSF/ANSI 49, by NSF International, a leading testing agency in North America and recognized globally.

Find out more at thermofisher.com/bsc

Australia +61 39757 4300 Austria +43 1 801 40 0 Belgium +32 9 272 54 82 **China** +800 810 5118, +400 650 5118 France +33 2 2803 2180

 $\textbf{Germany national toll free}\ 0800\ 1\ 536\ 376$

Germany international +49 6184 90 6000 India toll free 1800 22 8374 India +91 22 6716 2200

Italy +39 02 95059 552 Japan +81 3 5826 1616 Korea +82 2 2023 0600 $\textbf{Netherlands} + 31\ 76\ 579\ 55\ 55$ New Zealand +64 9 980 6700 Nordic/Baltic/CIS countries +358 10 329 2200

Russia +7 812 703 42 15, +7 495 739 76 41

Singapore +82 2 3420 8700

Spain/Portugal +34 93 223 09 18 **Switzerland** +41 44 454 12 12 UK/Ireland +44 870 609 9203 USA/Canada +1 866 984 3766

Other Asian countries +852 3107 7600 Countries not listed +49 6184 90 6000

