


# Sapphire™ Biomolecular Imager

THE NEXT GENERATION OF LASER-BASED IMAGING



# Breakthrough image capture and analysis


The Sapphire Biomolecular Imager is a next generation laser scanning system that provides you with exceptional data quality through extremely sensitive detection, ultra high resolution and broad linear dynamic range.



This system supports long and short wavelengths of near infrared fluorescence (NIR), red/green/blue (RGB) imaging, chemiluminescent imaging, phosphor imaging as well as optical densitometry (OD) of proteins in stained gels.\*

- Improved multiplex fluorescent detection
- Chemiluminescent imaging, surpassing film
- Higher sensitivity for lower limits of detection (femtograms)
- Broad linear dynamic range for accurate quantitation
- Ease of use with intuitive control software



\* Patent pending



## FLUORESCENCE IMAGING

### Improve Your Quantitation



- Stable fluorescent signals enable reliable quantitation of weak and strong bands
- Probe for multiple proteins on a single blot, and get more results from each sample



## CHEMILUMINESCENCE IMAGING

### Enhance your Detection


- Femtogram detection of protein
- The speed and sensitivity of traditional film with 4X the dynamic range



## PHOSPHOR IMAGING

### Improve Your Flexibility


- Scan storage phosphor screens for filmless autoradiography
- Wide dynamic range for sensitive detection



# FLUORESCENCE

Fluorescent probes are much more stable than traditional chemiluminescent methods since they do not rely on enzyme activity and substrate breakdown. This makes them ideal for quantitative analysis delivering consistent signal even months after the initial analysis.


Choose from multiple fluorescent dyes for sensitive results. Supporting a wide range of visible and NIR dyes, the Sapphire Biomolecular Imager enables low limits of detection.



Sample AzureSpectra 490 labeled purified BSA  
Imaging 488  
LOD 0.72 pg  
DR 6.3  
Linearity  $R^2=0.99$



Sample AzureSpectra 550 labeled purified transferrin  
Imaging 520  
LOD 2.9 pg  
DR 5.7  
Linearity  $R^2=0.992$



Sample AzureSpectra 650 labeled purified transferrin  
Imaging 658  
LOD 2.9 pg  
DR 5.7  
Linearity  $R^2=0.993$



Sample AzureSpectra 700 labeled purified transferrin  
Imaging 658  
LOD 1.45 pg  
DR 6.0  
Linearity  $R^2=0.999$




Sample AzureSpectra 800 labeled purified transferrin  
Imaging 785  
LOD 2.9 pg  
DR 5.7  
Linearity  $R^2=0.99$

## FLUORESCENCE: APPLICATIONS

### Quantitative Western Blots


Fluorescent Westerns enable quantification of post-translation modifications without stripping and reprobing. Spectrally separate the channels to resolve protein changes, and use a loading control to normalize for lane to lane loading variation.



Multiplex fluorescent detection of STAT1 (Cy3), Phospho-STAT1 (IR700), and GAPDH (IR800) in HeLa lysates. Composite image shows simultaneous detection of three proteins.

### See multiple proteins on the same blot


With multiplex fluorescence, you can detect up to 4 targets simultaneously, with a high degree of sensitivity in each channel. Reduce errors from blot to blot variations by probing for multiple targets on the same blot, without stripping and reprobing. With the Sapphire Biomolecular Imager, up to 4 different fluorescent probes can be used and imaged at once.



## FLUORESCENCE: APPLICATIONS


### High Resolution Gel Imaging

2D Electrophoresis and 2D-DIGE are used to accurately **quantitate subtle changes in protein expression**. High sensitivity, broad dynamic range, and resolution down to 10 microns make the Sapphire Biomolecular Imager suited for 2D electrophoresis applications, enabling you to detect changes in complex samples.




Untreated HeLa lysate was labeled with Cy3. Treated HeLa lysate was labeled with Cy5. Samples were simultaneously separated using IEF in the first dimension and SDS-PAGE in the second dimension.

**Laser-based imaging delivers sharp and crisp DNA and protein gel images.** Common DNA dyes such as EtBr and SYBR Green allow you to detect your DNA size or concentration. Visible protein stains, such as Coomassie Blue, and fluorescent protein stains, like Sypro Ruby, are also compatible with the Sapphire Biomolecular Imager.



Representative gels showing 3kb of DNA ladder.




Purified transferrin was serially diluted 1:1 from 5000 pg to 5 pg and separated using SDS-PAGE in a 4-15% Tris-Glycine gel. After separation, the gel was stained using Coomassie Blue and imaged using the 700 channel setting of the Sapphire Biomolecular Imager.

## FLUORESCENCE: APPLICATIONS

### Cell-Based ELISAs


Detect proteins in fixed cultured cells using target-specific primary antibodies and fluorescently conjugated secondary antibodies. Quickly and accurately measure relative protein expression in multiple samples, targeting several proteins with spectrally distinct dyes.



HeLa cells were serially diluted and seeded into a 96-well plate, cultured, fixed and permeabilized, then probed for Tubulin with AzureSpectra 550 (520 channel, green), beta-Actin with AzureSpectra 800 (785 channel, blue) and RedDot™1 Nuclear Stain as a normalization control (785 channel, red).

### Slide-Based Imaging

High resolution imaging, down to 10 microns, enables visualization of small features. Image fluorescent DNA and protein arrays, or tissue slides on the Sapphire Biomolecular Imager.



Full Moon BioSystems Scanner Calibration Slide scanned at 10  $\mu\text{m}$  resolution using the 520 nm and 658 nm channels of the Sapphire Biomolecular Imager.

# CHEMILUMINESCENCE



The sensitivity of chemiluminescent detection provides fast detection of low abundance proteins. Chemiluminescence can be more sensitive than other detection techniques due to the amplification of signal by enzyme activity. Additionally, most samples have no “background” signal, and are ideal for researchers who want the cleanest images.




# CHEMILUMINESCENCE: APPLICATIONS

## Western Blotting

Digital imaging enables a wide dynamic range for chemiluminescent imaging. Traditional film saturates quickly, making it hard to distinguish weak and strong bands simultaneously.



Purified transferrin serially diluted and detected with chemiluminescent Western using Azure Radiance Substrate. Exceptional signal over background ratios allow for detection of small quantities of protein with the Sapphire Biomolecular Imager compared to film.




Simultaneous imaging of chemiluminescent samples and colored molecular weight markers. Images are captured serially, and then overlaid within the capture software.


# PHOSPHOR IMAGING

To detect radioactive signals using phosphor imaging, samples containing radiolabeled samples are exposed to a storage phosphor screen. Light is emitted from the screen in proportion to the amount of radioactivity in the sample upon laser-induced stimulation.

## Filmless autoradiography



Sample <sup>14</sup>C autoradiographic standard  
Imaging Ex 658 / Em 488  
LOD 0.036  $\mu$ Ci/g  
DR 5.4 orders of magnitude  
Linearity R<sup>2</sup>=0.99




American Radiolabeled Chemicals Carbon-14 standard exposed to Storage Phosphor Screen for three hours then imaged on the Sapphire Biomolecular Imager.

# IMAGE CAPTURE AND ANALYSIS

## Workflow software for gels and blots

Application driven software makes setting up imaging protocols **simple**. Sapphire Capture software is designed to walk you through the steps to set up the perfect imaging parameters.




AzureSpot Analysis software designed for quantitation. View and analyze single and multichannel images.




# TECHNOLOGY DRIVEN PERFORMANCE

Discover more with a laser-based imager. The quality of excitation light affects image quality and your ability to detect subtle changes. The Sapphire Biomolecular Imager uses up to 4 solid state lasers as excitation sources.




# 1 PHOTON MULTIPLIER TUBE (PMT)

Highly sensitive photodetector




# 2 AVALANCHE PHOTODIODE (APD)

High quantum efficiency in the NIR



# 3 CHARGE-COUPLED DEVICE (CCD)

Samples the entire spectrum simultaneously, and allows long integration times



The Sapphire Biomolecular Imager's patent pending design uses 3 different detection modes to deliver ultimate application flexibility. Using PMTs, APDs, and CCD, the Sapphire is a unique system designed for visible fluorescence, NIR fluorescence, and chemiluminescent imaging.



## Choose from three models

### SAPPHIRE NIR

---

NIR Fluorescence	✓
Chemi Imaging	optional
Phosphor Imaging	optional

### SAPPHIRE RGB

---

RGB Fluorescence	✓
Chemi Imaging	optional
Phosphor Imaging	optional

### SAPPHIRE RGBNIR

---

RGB Fluorescence	✓
NIR Fluorescence	✓
Chemi Imaging	optional
Phosphor Imaging	optional

## Specifications

Part number

IS1000

IS1001

IS1002

Laser excitation wavelengths

658, 785

488, 520, 658

488, 520, 658, 785

Bit depth

16 bit

16 bit

16 bit

Scanning area

25 cm X 25 cm

25 cm X 25 cm

25 cm X 25 cm

Scanning mode

Simultaneous

Simultaneous

Simultaneous

Pixel size

10 micron – 1000 micron

10 micron – 1000 micron

10 micron – 1000 micron

Detectors

Avalanche Photodiode

Photomultiplier Tube, Avalanche Photodiode

Photomultiplier Tube, Avalanche Photodiode

Scanning speed

50 cm/s

50 cm/s

50 cm/s

Filters

Red 721BP13, IR 832BP37

Blue 515BP15, Green 565BP12,  
Red 721BP13

Blue 515BP15, Green 565BP12,  
Red 721BP13, IR 832BP37

Dimensions (W X H X D)

75 cm X 45 cm X 70 cm

75 cm X 45 cm X 70 cm

75 cm X 45 cm X 70 cm

Weight

140 lbs

140 lbs

140 lbs

Power requirements

100 – 240 VAC ± 10%, 50/60 Hz

100 – 240 VAC ± 10%, 50/60 Hz

100 – 240 VAC ± 10%, 50/60 Hz

Computer

Windows 10 laptop

Windows 10 laptop

Windows 10 laptop

**Chemiluminescent Upgrade**

**Optional**

**Optional**

**Optional**

Imaging area, chemiluminescent

16 cm X 13 cm

16 cm X 13 cm

16 cm X 13 cm

Bit depth

16 bit

16 bit

16 bit

Resolution

2688 X 2200

2688 X 2200

2688 X 2200

Imaging time

0-60 minutes

0-60 minutes

0-60 minutes

Visible imaging

Yes

Yes

Yes

**Phosphor Imaging Module**

**Optional**

**Optional**

**Optional**

Sample types

Storage phosphor screen (imaging plate)

Storage phosphor screen (imaging plate)

Storage phosphor screen (imaging plate)



6761 Sierra Court, Suite F • Dublin, CA 94568 USA

Phone: 925.307.7127 • Fax: 925.905.1816

info@azurebiosystems.com • www.azurebiosystems.com

Copyright © 2017 Azure Biosystems. All rights reserved. The Azure Biosystems logo, Azure Biosystems™, and Sapphire™ are trademarks of the Company. Alexa Fluor® and QDots® are a trademark of Molecular Probes, Inc. Cy3®, Cy5® and Cy2® are registered trademarks of Amersham Biosciences. DyLight® is a trademark of Thermo Fisher Scientific Inc. RedDot™ is a trademark of Biotium. All other trademarks, service marks and trade names appearing in this brochure are the property of their respective owners.