

3D TableTrix® Microcarriers

Porous & Dissolvable Pharmaceutical-grade Microcarriers
for Adherent Cell Mass Manufacturing

Architect for Cells

— Expert in 3D manufacturing of high quality cells

Cell Mass Manufacturing

Who Needs It & How It is Done

Biotech Sectors in Need of Cell Mass Manufacturing



Vaccine



Antibody



Oncolytic Virus

Biologics



Cell Therapy



Tissue Engineering



Gene Therapy

**Regenerative
Medicine**



Extracellular
Vesicles



Secretomes

Cell-Free Therapy



Cultured Meat




More Potential
Applications...

**Future Food
Tech**

Options for Cell Manufacturing

2D Culture




Common cell types:

Most adherent cells

Industrial applications:

- Stem cell therapy
- Diploid cells for vaccine production
- Virus packaging with HEK293T
- ... others that current suspension or microcarriers cannot support

Suspension Culture




Common cell types:

Adapted suspension mammalian cells (CHO, 293F, MDCK)
Insect cells
Immune cells (T cells, NK, etc)

Industrial applications:

Biologics/Vaccines & proteins
Gene therapy/Virus packaging
Immune cell therapy

Microcarriers & Fibers




Common cell types:

Adherent mammalian cells
(293, adherent CHO, VERO, MARC145, immortalized hepatic cells, etc.)

Industrial applications:

Biologics/Vaccines & proteins
Gene therapy/Virus packaging
Exosomes

Limitations






3D TableTrix® Microcarriers

Dissolvable porous microcarriers

What is Microcarrier

- A micron-size support matrix that allows for the growth of adherent cells in bioreactors
- Usually exist as 100–300 μm polymeric beads



3D TableTrix® Microcarriers: Unique Features



Dissolvable

>99% recovery of cells



High Porosity

Up to 9000 cm²/g
Yield 10⁸~10⁹ cells



Pharmaceutical Grade





Produced under GMP
Registered with U.S. DMF



Ready-to-use

Sterile, weight-defined*

*Applicable to some products only






Dissolvable
>99% recovery of cells


Unique Features: Dissolvable

Complete cell recovery from microcarriers
by dissolving with 3D FloTrix® Digest



Performed at 37°C, pH 7.0-7.4

Cells remain viable during process



LIVE/DEAD


Unique Features: Dissolvable

Low residual content in end products

— remove simply by centrifuge



Accredited by National Institutes for Food & Drug Control (China)



No toxicity was observed in animals injected with 1X, 10X or 100X concentrations of residual content thresholds



Dissolvable: Why is It Important?

- Cells are THE **“living” end-product** of cell therapy
 - Downstream processing methods for live cells are limited
 - Extensive purification cannot be performed
 - Need to keep residuals from starting materials used in the manufacturing process low to ensure safety
- Good **quality seed cells** ensure efficient biologic production
 - Large quantity of cells needed for the final stage of the upstream process (i.e. virus production, protein production) in the largest bioreactor (up to 2000 L bioreactors)
 - Harsh detachment methods used in conventional undegradable microcarriers are not suitable for delicate cells, i.e. diploid cells
 - Biomimetic environment to trigger better virus production



Ready-to-use
Sterile, weight-defined*

Unique Feature: Dispersible Tablets

- Ready to use



- Pre-sterilized
- Closed system packaging for bioreactors (* powder only)

- Designed for lab-scale process

- Defined weight per tablet (20 mg)
- No hassle of weighing powder
- Static culture in 6-well plates (one tablet/well)
- Dynamic culture in spinner flasks or with microSPIN* system

- Dispersible

- **NO** adhesive or additional chemical substances
- **NO** disintegrants
- Tablet disperses into individual microcarriers when in contact with liquid



* Bioreactor designed with 6-well plates with built-in impellers, available from CytoNiche

Unique Feature: Safe for cellular drugs



Pharmaceutical Grade
Registered with U.S. DMF



Approved by FDA

DMF # 35481

One and only microcarriers listed by FDA under DMF filings



Pharmaceutical Excipients Recognized by the Center for Drug Evaluation (China)

Registration No.

#F20200000496

#F20210000003



Quality and Safety Report

- ✓ Biocompatibility
- ✓ Stability
- ✓ Cytotoxicity
- ✓ Acute Toxicity (*in vivo*)
- ✓ Allergic reaction
- ✓ Immunogenicity
- ✓ Hemolysis test
- ✓ Pyrogenic reaction
- ✓ ...

List of Drug Master Files (DMFs)

The list of DMFs, which is updated quarterly, contains DMFs RECEIVED by March 31, 2021, for which acknowledgment letters were sent before April 8, 2020.

The list is current through DMF 035821. Changes to the DMF activity status (A=active; I=inactive), DMF type, holder name and subject (title) made since the last update of December 31, 2020, are included.

• 1Q2021 - All - Excel (XLS)

DMF#	STATUS	TYPE	SUBMIT DATE	HOLDER	SUBJECT
35481	A	IV	3/1/2021	CYTONECHE	GELATIN MICROCARRIER TABLET FOR CELLS, STERILE, EXCIPIENT

国家药品监督管理局药品审评中心
CENTER FOR DRUG EVALUATION, NMPA

原料药、药用辅料和药包材登记信息公开

登记号	品种名称	企业名称	企业地址	产品来源	包装规格	规格	更新日期	与辅料及药包材关系	备注
F2020000003	微载体	天津康希诺生物股份有限公司	天津滨海新区经济技术开发区第十大街104号10407室	国内	100片/瓶	100片/瓶	2020-02-04	1	
F20200000496	微载体	天津康希诺生物股份有限公司	天津滨海新区经济技术开发区第十大街104号10407室	国内	100片/瓶	100片/瓶	2020-12-07	1	

总结报告


明胶微载体溶血试验

试验名称: 微载体溶血试验
 批号: KY2010104
 试验负责人: 刘慧
 试验日期: 2020.06.25-2020.06.30

试验结果: 微载体在生理盐水、蒸馏水、0.9%氯化钠溶液中均无溶血现象, 符合药典要求。

3D TableTrix® Microcarriers: Comparison

Conventional Microcarriers



- ❑ Non-degradable
 - ❑ Harsh harvesting condition
 - ❑ Extra separation protocol
- ❑ Pre-treatment needed
 - ❑ Powder in bulk only
 - ❑ May need bead hydration
- ❑ Solid sphere or macroporous
- ❑ Synthetic polymers




3D TableTrix® Microcarriers



- ✓ Fully dissolvable
 - ✓ Gentle harvesting
 - ✓ No extra separation
- ✓ Ready-to-use*
 - ✓ Weight-defined tablets
 - ✓ Dispersible
- ✓ Macroporous
- ✓ Biomimetic protein

*Choose from:


- Ready-to-use weight-defined dispersible tablet for lab scale
- Powder in bulk for mass production (closed system packaging)



3D TableTrix® Microcarriers: Comparison


Product / Feature	Cytodex-1	Cytodex-3	SoloHill	Cytopore1	Cytopore2	3D TableTrix W01	3D TableTrix G02	3D TableTrix V01
Material	Dextran	Dextran	Polystyrene	Cellulose	Cellulose	Denatured collagen	Denatured collagen	Denatured collagen
Structure	Solid Sphere	Solid Sphere	Solid Sphere	Porous Sphere	Porous Sphere	Porous Sphere	Porous Sphere	Porous Sphere
Size (µm)	147-248	133-215	90-150/125-212	200-270	200-270	100-200	120-180	80-200
Pore (µm)	/	/	/	30	30	30	30	30
Surface Area (cm²/g)	4400	2700	360/480	~11000	~11000	>6000	>6000	>6000
Surface Modification	DEAE positive charge	Gelatin	RGD/Positive charge/Gelatin	DEAE positive charge	DEAE high positive charge	N/A	N/A	N/A
Packaging	Powder	Powder	Powder	Powder	Powder	Sterile Tablet/Powder	Sterile Tablet/Powder	Powder
Harvest Efficiency	Poor	Poor	Poor	Not possible	Not possible	High	High	High
Degradability	NO	NO	NO	NO	NO	YES	YES	YES
Application Field	Biologics	Biologics	Biologics	Biologics	Biologics	Cell & Gene Therapy Biologics	Cell & Gene Therapy Biologics	Biologics

3D TableTrix® Microcarriers: Applications

Product / Feature	3D TableTrix W01	3D TableTrix G02	3D TableTrix V01
Material	Denatured collagen	Denatured collagen	Denatured collagen
Size (µm)	100-200	120-180	80-200
Packaging	Tablet/Powder	Tablet/Powder	Powder
Pre-sterilized	Yes	Yes	No
Suitable cell types	 Mesenchymal stem cells	VERO MDCK 2BS/MRC-5 Fish primary cells	HEK293T Adherent CHO
Suitable application	Working cell bank expansion Final cell production Exosome production	Seed cell expansion Human/Animal vaccine production Oncolytic virus production	Lenti-virus production Adeno-associated virus production Protein expression

Cell Mass Manufacturing Process with 3D TableTrix®

Use of 3D TableTrix® microcarriers for adherent cell manufacturing



Cell Types We Have Worked With

Cell Therapy

hUC-MSC	Placenta MSC
Adipose-derived MSC	Amniotic membrane MSC
Dental Pulp MSC	ESC-derived MSC
Bone Marrow MSC	iPSC-derived MSC

Vaccines

MRC-5	Vero
2BS	PK-15
Diploid cells	MARC-145
MDCK	Fish cells

Virus Packaging

HEK293T cell line

Artificial Liver

Immortalized Human
Hepatocytes

Other Cell Types



hAEC (Human amniotic epithelial stem cells)	HDF Human Dermal Fibroblasts
UPC human kidney adult stem cells	3T3 Cell Lines
Nasal mucosal stem cells	Neurons
HUVEC	

Cultured Meat

Porcine Muscle Stem Cells
Porcine Adipose Stem Cells

Performance Data: MSCs with W01 microcarriers


A




Serum-free medium, 4-5 days culture

Data from Journal of Tissue Engineering and Regenerative Medicine, 2022, 16(10), 934-944.

hUC-MSC Expansion



3D FloTriX®
miniSPIN



Phase 1

3D FloTriX®
vivaSPIN (5L)



Phase 2

3D FloTriX®
vivaSPIN (15L)



Phase 3


3D FloTriX®
vivaPREP PLUS



Wash & Concentrate

Check MSC deck for more details on cell quality & performance comparison with other microcarriers



Performance Data: V01 microcarriers






Performance of V01 Microcarrier

10^7 cells/mL is achievable
Performance validated in 300 L bioreactor
Bead-to-bead transfer for passaging

Check Vaccine & Virus deck for more details on more vaccine virus titer & performance comparison with other microcarriers & 2D



Performance Data: HEK293T on G02 microcarriers



Performance of G02 Microcarrier

10⁷ cells/mL of HEK293T is achievable
Density comparable to suspension culture
Bead-to-bead transfer for passaging

Check HEK293T deck for more details on more plasmid transfection efficiency, lenti-virus & adeno-associated virus titer, & cell culture performance comparison with other microcarriers

Performance Data: various cell types

Cells	Culture Medium	Culture time (days)	Expansion Factor	Cell Viability
MSCs from various tissue sources	Serum-free	4-6	8-14	≥95.7%
Porcine muscle stem cells	Serum-free	7	16.52	95.1%
Human Hepatic Cells	Serum-free	7	2.8	90.3%
HUVEC	With serum	10	35	95.5%
3T3	With serum	6	24.55	92.6%
293T	With serum	4	9.3	94.0%
PK15	With serum	3	5.4	90.7%
Vero	With serum	4	8.2	90.4%
MRC-5	With serum	4	5.0	97.47%

Holistic Solution for Cell Mass Production

CytoNiche offers various products & extensive technical support to accelerate your cell manufacturing process

Reagents and Consumables




3D TableTrix®
Microcarriers



3D FloTrix®
Digest for microcarriers
3D FloTrix®
MSC SF Medium

Bioreactors



3D FloTrix®
micro/miniSPIN
Bioreactors

3D FloTrix®
vivaSPIN
Stir-tank
Bioreactors



Single-use
3D FloTrix®
vivaROCK
Bioreactors

Downstream Processing



3D FloTrix®
vivaPREP/vP PLUS
cell processing system



3D FloTrix®
vivaEXO
Exosome Collection System

Holistic Solution for Cell Mass Production


New products are coming Sep 2023 to bring cell mass production at GMP-standard a step further.

Microcarriers



**Animal
Component Free
Dissolvable
Microcarriers**

Bioreactors



3D FloTrix® miniSPIN FLEX

Independent control for each channel
Ultra slim to fit 2 sets into an incubator
with disposable spinner flasks

Downstream Processing



3D FloTrix® vivaPACK
cell formulation
& filling system
(cryobag)/(cryovial)

Publications

- **Dispersible and dissolvable porous microcarrier tablets enable efficient large-scale human mesenchymal stem cell expansion (2020).** *Tissue Engineering Part C: Methods*, 26(5), 263-275.
- **GMP-grade microcarrier and automated closed industrial scale cell production platform for culture of MSCs (2022).** *Journal of Tissue Engineering and Regenerative Medicine*, 16(10), 934-944.
- Engineering 3D functional tissue constructs using self-assembling cell-laden microniches (2020). *Acta biomaterialia*, 114, 170-182.
- Intra-articular injection of cell-laden 3D microcryogels empower low-dose cell therapy for osteoarthritis in a rat model (2020). *Cell transplantation*, 29, 0963689720932142.
- A low dose cell therapy system for treating osteoarthritis: In vivo study and in vitro mechanistic investigations (2022). *Bioactive Materials*, 7, 478-490.
- Exendin-4 gene modification and microscaffold encapsulation promote self-persistence and antidiabetic activity of MSCs (2021). *Science Advances*, 7(27), eabi4379.
- Cell-subpopulation alteration and FGF7 activation regulate the function of tendon stem/progenitor cells in 3D microenvironment revealed by single-cell analysis (2022). *Biomaterials*, 280, 121238.
- Efficient endothelial and smooth muscle cell differentiation from human pluripotent stem cells through a simplified insulin-free culture system (2021). *Biomaterials*, 271, 120713.
- Engineered meatballs via scalable skeletal muscle cell expansion and modular micro-tissue assembly using porous gelatin micro-carriers (2022). *Biomaterials*, 287, 121615.
- Exosomes derived from three-dimensional cultured human umbilical cord mesenchymal stem cells ameliorate pulmonary fibrosis in a mouse silicosis model (2020). *Stem Cell Research & Therapy*, 11(1), 1-12.
- Exosomal let-7i-5p from three-dimensional cultured human umbilical cord mesenchymal stem cells inhibits fibroblast activation in silicosis through targeting TGFBR1 (2022). *Ecotoxicology and Environmental Safety*, 233, 113302.

Empowering cell and gene therapy with intelligent 3D cell mass manufacturing technology to benefit more patients

Architect for Cells

— Expert in 3D manufacturing of high quality cells



YAN XIAOJUN

Co-founder & CTO

WhatsApp: 65-8698-0015

email: yanxiaojun@cytoniche.com

CytoNiche SG

Showroom: To be coming Sep 2023

Launching a new era for industrialising cell manufacturing

Headquarter

B1-208 Yiyuan Culture & Innovation Base, Haidian, Beijing, P.R.China

 86-010-6286 8280

 marketing@cytoniche.com

 en.cytoniche.com

Collaboration & Support Center

Blk 4, 260 Maoyuan Road, Fengxian District, Shanghai, P.R.China

GMP Production Center

402, Plant 6, 1 Xinxing Road, Wuqing District, Tianjin, P.R.China