

# Introduction of 3D FlorTrix<sup>®</sup> vivaEXO

Architect for Cells

— Expert in 3D manufacturing of high quality cells

---

# CONTENTS

Solidarity

High quality and efficiency

Innovation and creation

## 1. CytoNiche's Core Technology

## 2. 3D FloTrix® vivaEXO

2-1 Introduction

2-2 Performance data

2-3 Summary


# CytoNiche's Core Technology

Cell culture technology based on 3D TableTrix<sup>®</sup> dissolvable microcarrier  
and 3D FloTrix<sup>®</sup> process.



# Core Technology: 3D FloTrix® Cell Culture Technology

Using 3D TableTrix® dissolvable microcarriers for adherent cell manufacturing




# 3D FlorTrix® vivaEXO

For Supernatant Clarification, Purification & Concentration



# 3D FloTrix® vivaEXO - Exosome Harvesting System



Dimensions (W\*L\*H) 1500×800×1000mm

Weight 40kg



Large Capacity



Automated & Continuous  
& Semi Closed




High Efficiency



GMP Compliances




Principle: Multi-stage Filtration





## 2-1 3D FloTrix® vivaEXO - Working principle



Process steps:

① Activation (optional)

--- Sample processing ---

② Clarification and filtration: CLR1 and CLR2

-R1 (Intermediate product bag)

③ Concentration filtration

-The permeate enters R2 (waste bag), retentate returns to R1


④ Sterile filtration → R3 (product bag)

⑤ Cleaning (optional)

## 2-1 3D FloTrix® vivaEXO

### 1. Easy to use

- Clearly labelled tubes
- Conveniently installed
- Automated



Standard consumable kit

### 2. Closed processing

- Fully closed after connecting with luer connectors
- Includes final sterile filtration



Consumable kit installation



**R1, R2, R3 fluid storage bags:**  
**endotoxin <0.125EU/mL**



## 2-1 3D FloTrix® vivaEXO

### 3、 Scale processing

- Process **10L** of starting solution per batch(**1-3 hours**).
- 10L starting solution → 300~500mL final product  
(**concentrated ~20 times**).

### 4、 Automated processes controlled by mass and pressure sensors

- Mass sensor:
  - Associated with peristaltic pump 2 to start concentration filtration
  - Real-time monitoring of filtration progress
- Disposable pressure sensor:
  - Associated with the peristaltic pump 2 to monitor transmembrane inlet pressure
  - Prevents high pressure from damaging membrane & breaking delicate extracellular vesicles



## 2-1 3D FloTrix® vivaEXO

### 5、 Scalable

- Easily realize large-scale processing of cell supernatant harvest (10L-30L, only by changing CLR module every 10L)
- Concentrate and diafiltrate to improve product purity



### 6、 GMP Compliant

- User login and multi-level permission management
- Exportable records
- Data integrity and traceability




## 2-2 3D FloTrix® vivaEXO - Performance data

### Application ①



30000X  
200nm






- ✓ Particle size: 50-200nm
- ✓ Exosome protein markers: clear bands of CD9, CD63, and TSG101, negative for Calnexin
- ✓ Obvious concave structure observed under TEM


Sample source : Customer case using UC-MSC supernatant harvested from 2D culture

## 2-2 3D FloTrix® vivaEXO - Performance data

### Application ②



与细胞和细胞外基质共存  
H-7650 30 kV 200000 x



与细胞和细胞外基质共存  
H-7650 30 kV 200000 x

- ✓ Particle size: 50-150nm
- ✓ Exosome protein markers: clear bands of TSG101, CD81, and CD63, negative for Calnexin
- ✓ Obvious concave structure observed under TEM

Sample source : UC-MSC supernatant harvested from 3D FloTrix® process




## 2-2 Application of 3D FloTrix® Technology - Stem cell exosome preparation

### Cell preparation



Day 0~3 Cell cultivation in 2D

### 3D TableTrix® microcarrier 3D FloTrix® miniSPIN



Day 3~7 Cell cultivation in 3D

### 3D TableTrix® microcarrier 3D FloTrix® vivaSPIN



Day 7~11 Supernatant harvest

### 3D FloTrix® vivaEXO Exosome Harvesting System



Day 11 Exosome concentration  
& purification

### Process steps:

Phase I: 2D recovery of cells (72h)

Phase II: miniSPIN system for 3D culture of seed cells (96h)

Phase III: vivaSPIN-10L for 3D culture cells and produce exosomes on a large scale (96h)

Phase IV: Use vivaEXO to purify and concentrate exosomes (1~3h)



## 2-3 Summary

### Advantages of harvesting exosomes using 3D FloTrix® technology.

#### 3D biomimetic culture

Simulate the living environment of cells  
Promote cells to secrete high-quality exosomes

#### Reliable quality

Collect exosomes and harvest cells at the same time  
Meet quality inspection requirements.

10g microcarriers + 10L medium  
+ cells

CytoNiche  
Architects for Cells

3D FloTrix®  
Technology

#### Short cultivation period

A single batch is only 4 to 5 days

#### High Yield

Harvest  $10^{13} \sim 10^{14}$  exosomes  
single batch

#### Support customization

Support 3D FloTrix® customized  
process to harvest exosomes

4~5 days

$2 \sim 7 \times 10^{13}$  exosomes

\*amount of exosomes collected highly depends on medium used, data shown is generated from RMZ112

## Product Bundles

Category No.	Product Name	Qty	Spec	Notes
FTVE10	3D FloTrix® vivaEXO Exosome Harvesting System	1 unit	Controller 1unit, GL45 Bottle Cap with 3-ports and tubes* 3sets	NA
R010-CLR-01	3D FloTrix® vivaEXO Exosome CLR Processing Kit	1 cas	2set/cas	NA
R011-PUR-02	3D FloTrix® vivaEXO Exosome PUR02 Processing Kit	1 cas	1set of PUR /cas, small MWCO	Single use
R011-PUR-03	3D FloTrix® vivaEXO Exosome PUR03 Processing Kit	1 cas	1set of PUR /cas, medium MWCO	3 times use within a month

# THANKS



Beijing Cytoniche Biotech.Co.,Ltd.  
Room.208, Block B1, Yiyuan Cultural Base District B, Haidian District, Beijing.  
[400-012-6688](tel:400-012-6688) | [sales@cytoniche.com](mailto:sales@cytoniche.com) | [www.cytoniche.com](http://www.cytoniche.com)

Shanghai CytoNiche Biotech Co., Ltd.  
15F, Building 4, No.260 Maoyuan Road, Fengxian District, Shanghai  
[\(86\)021-61905266](tel:86021-61905266)