

Ensure **ADDITIONAL**
LAB SAFETY during
SARS-Cov-2 inactivation
with Heratherm APS
incubator Auto-
Decontamination Function



140 °C push button decontamination



Building on our established CO₂ incubator decontamination technology, introducing the first microbiological incubators with an independently certified 140 °C decontamination routine.

140 °C decontamination

- At 140 °C contaminating microorganisms are reduced to a minimum, comparable to sterilization, within a six hour cycle
- No need for separate autoclaving of interior fittings
- Certified by an accredited microbiological institute

Superior sample protection

- Automatic over-and under temperature alarm

Intelligent design for improved results

- Lockable incubator door for restricted access
- Audible alarm if the door is left open accidentally



[Click here to know more](#)

Find out more at: **thermofisher.com**

General Enquiries: **Laboratory.Products.SEA@thermofisher.com**

Singapore Office:

Thermo Fisher Scientific
33 Marsiling Industrial Estate Road 3
#07-06 Singapore 739256

Taiwan Office:

台灣賽默飛世爾科技股份有限公司
11493 台北市內湖區洲子街 85
號 6 樓之1

ThermoFisher
SCIENTIFIC