

# MOBILE LABORATORIES

Covid-19 Mobile Response






**Keeping Us Safe.**









## Biological Mobile Laboratory

Typical BSL 2 Lab



Front View

# General Mobile Laboratory
















Front View





### Colors Available

	
White	Purple White
	
Cream	Light Grey
	
Dark Grey	Black
	
Traffic Red	Orange
	
Signal White	Pastel Green
	
Emerald Green	Light Blue
	
	Dark Blue





**Stainless Steel Fume Hoods  
with Built-in Scubbers for Safety**



**Thermo Scientific**  
Class II Type A2



**Baker IsoGuard**  
Class III Glovebox

# Specifications

## Design Flexibility

- 20ft / 40ft Container / Multiple configurations
- Biosafety Level 2 / 2+
- 10mm thick lab grade phenolic walls - easily disinfected
- Vinyl sheet flooring with coved corners
- LED lights, delivering up to 800 lux
- Lab grade AdvancelabPro<sup>®</sup> Casework,
- quality steel construction
- Lab grade phenolic resin table with backsplash
- Lab sink and tap, lever type
- Built-in air conditioner system
- Electrical distributor board with MCCB for hook up at site (30A DB Single Phase)



# Options



Auto Sash /  
Self-closing Controller



Stairs / Ramps



Custom Paint Logo



Pass Through Box



Portable Satellite Data Antenna



Field Clean Booth



Generator



Negative Air  
Pressure Package

# Our Clients



And many more...

# Enquiry Form

Once we have received the filled form, we will contact you shortly.\*

Full Name:

Salutation: Dr / Mr / Mdm / Ms

Email:

Date:

Company:

Contact Number:

Expected Delivery Date :

Requirement / Usage:

Please describe your specifications:

\*Please kindly send your enquiry to [info@advancelab.com.sg](mailto:info@advancelab.com.sg)



# Warranty

**Products Manufactured by Advancelab:** Advancelab(S) Pte.Ltd., warrants products that it manufactures to be free from defects for a period of 12 months for parts, commencing from the date of shipment. Advancelab's sole responsibility is to repair or replace, at its option, any part of the product that proves defective or malfunctioning during this time limit. This warranty is void if the equipment is abused or modified by the customer, is operated outside Advancelab's operating instructions or specifications, or is used in any application other than that for which it is specified. This warranty does not include routine maintenance or service procedures, breakage, shipping damage, nor damage from misuse, intentional or unintentional abuse, neglect, natural disasters, or acts of God.

**Freight Shortage or Damage:** Upon receipt of any equipment from Advancelab, customer shall immediately unpack and inspect for damage or shortage. The customer shall not accept a damaged package or a short shipment until the carrier makes a "damage or shortage" notation on both the carrier's and customer's copy of the freight bill or delivery receipt. Service title passes when the shipment is loaded, so customer is responsible for filing and collecting a freight claim. Any replacement products must be ordered and paid for separately.

Generally, customers can improve the chance of collecting on a freight claim by following these procedures:

1. Formally requesting that the carrier inspect the shipment immediately upon suspecting damage or shortage to verify condition.
2. Notifying the carrier upon discovery of concealed damage and requesting an inspection within 15 days of receipt, both in person or phone and following up via mail.
3. Keeping the shipment as intact as possible, including retaining original packaging materials and keeping the product as close to the original receiving location as possible.
4. Holding salvage for disposition by the carrier.

**All Claims:** Advancelab(S) Pte.Ltd., expressly disclaims all other warranties, expressed or implied or implied by statute, including the warranties of merchantability or fitness for intended use. Advancelab is not responsible for consequential or incidental damages arising out of the purchase or use of the products supplied by Advancelab. Advancelab is not liable for damage to facilities, other equipment, products, property or personnel of others, or of their agents, suppliers, or affiliated parties, which is caused or alleged to have been caused by products supplied by Advancelab. In any event or series of events, Advancelab's total liability for any and all damages whatsoever is limited to the lesser of the actual damages or the original invoice cost of the items alleged to have caused the damage. The customer's sole and exclusive remedy for any cause of action whatsoever is repair or replacement of the non-conforming products or refund of the actual purchase price, at the sole option of Advancelab. All claims must be made in writing within 90 days of the date the product was shipped. Any claims not made within this time limit shall be deemed waived by the customer. Advancelab is not responsible for any additional costs of repair caused by poor packaging or in-shipment damage during return.

**Warranty Returns:** All warranty returns must be authorized in advance by Advancelab and approved by writing. Unless approved in advance for good reason, all returns must be in original condition, including all manuals, and must be packaged in original packaging materials. All returned goods are to be shipped to Advancelab, freight prepaid at customer's expense.

## Contact Details

### Singapore

Advancelab (S) Pte Ltd

253 Kaki Bukit Ave 1,  
Singapore 416061

T: +65 6448 8255  
F: +65 6448 9833  
E: [info@advancelab.com.sg](mailto:info@advancelab.com.sg)

### Malaysia

Advancelab Sdn Bhd

Johor Factory: No. 3388, Jalan  
Pekeliling Tanjung 27/2, Kawasan  
Perindustrian Indahpura, 81000  
Kulaijaya, Johor, Malaysia.

T: +607 660 8877  
F: +607 660 8866  
E: [info-my@advancelab-global.com](mailto:info-my@advancelab-global.com)

### Myanmar

Advancelab Scientific &  
Engineering Co., Ltd

No.(976/Ka), Than Thu Mar Road,  
5th Ward, South Okkalapa Township,  
Yangon, Myanmar.

T: +95 (1) 572393 /  
+95 (9) 779753802  
F: +95 (1) 572393  
E: [info-mm@advancelab-global.com](mailto:info-mm@advancelab-global.com)

### Indonesia

PT Advancelab Saintifik

253 Kaki Bukit Ave 1,  
Singapore 416061

T: +65 6448 8255  
F: +65 6448 9833  
E: [info@advancelab.com.sg](mailto:info@advancelab.com.sg)

### Thailand

Prima Scientific Co., Ltd.

147 / 170 - 171 Baromrajchonnanee  
Road, Arunamarin Bangkoknoi,  
Bangkok 10700, Thailand.

T: +66 2 884 9480  
F: +66 2 884 6441  
E: [primasci@primasci.com](mailto:primasci@primasci.com)

### Website:

[www.advancelabglobal.com](http://www.advancelabglobal.com)

