

Thermo Scientific Locator PLUS

Liquid N₂ Cryogenic Storage Systems

Maximum storage capacity with minimum liquid nitrogen consumption

Order Information

Product Name	Locator 4 PLUS	Locator 6 PLUS	Locator 8 PLUS
Model No. (without level monitor)	CY50935-70	CY50985-70	CY50945-70
Model No. (With level monitor 100-240V)	CY509108	CY509109	CY509111
Storage Type	Manual Fill / Rack & Box	Manual Fill / Rack & Box	Manual Fill / Rack & Box
LN ₂ Capacity	121 Liters	184 Liters	121 Liters
Neck Diameter	8.5 inch / 21.5cm	8.5 inch / 21.5cm	6.0 inch / 15.2cm
Static Evaporation Rate	0.99 Liters / Day	0.99 Liters / Day	0.60 Liters / Day
Static Holding Time	122 Days	185 Days	201 Days
Number of 1.2 / 2mL vial racks	4	6	8
Number of 2mL vial boxes per rack	10	10	10
Number of 2mL vials per box	100	100	25
Vial Capacity 2mL	4,000	6,000	2,000
Vessel Exterior Dimension			
Diameter x Height	55.8 x 95.3cm	66.0 x 95.3cm	55.8 x 95.3cm
Optional Accessories			
Low Level Alarm (stand alone)	AY509X1	AY509X1-70	AY509X1
Wheeled Accessory Cart (5 inches high)	CN509X4-70	CN509X11-70	CN509X7-70
CryoClaw	AY509X18	AY509X18	AY509X18
LN2 Dip Measurer	180143	180143	180143
Customer installed level monitor 100-240V	1950338	1950338	1950338



Freeze cell using tubes from 1 to 2ml. (up to 18 tubes) using Mr.Frosty freezing container. The simple to use system is designed to achieve a rate of cooling very close to -1°C/minute, the optimal rate of cell preservation

Special Offer

Purchase Locator PLUS get 2 units of "Mr.Frosty" freezing container and free accessories below:

Wheeled Accessory Cart (5 inches high)
CryoClaw
LN2 Dip Measurer



Contact your account representative for more information.

บริษัท ปริมา ไซเอนติฟิค จำกัด 147/170-171 ถ.บรมราชชนนี แขวงจตุจักร เขตจตุจักร กรุงเทพฯ 10700

Tel / โทรศัพท์ 02-435-7434, 02-434-3771, 02-433-2081 Fax / โทรสาร 02-884-6441 Email / อีเมล primasci@primasci.com



Thermo
SCIENTIFIC