

Thermo Scientific Refrigerated Incubators Your samples, safe inside



From incubation and sample storage in life sciences to BOD (biochemical oxygen demand) testing in water treatment plants, the need for temperature precision is critical. The Thermo Scientific™ Refrigerated Incubators, designed to fit your budget, are quality products with features that provide precision and produce consistent results − bringing basics to a new level.



standard 30 to 37 °C incubation applications in warm environments.

Thermo Scientific Refrigerated Incubators

Thermo Scientific Refrigerated Incubators offer an array of features and add-on accessories designed to provide users with ease-of-use and protection for valuable samples.









Model	RI-150 (Table Top)	RI-250 (Floor Standing)
Order Number	RI-150CN	RI-250CN
Chamber Volume	150 L / 5.3 cu ft	250 L / 8.8 cu ft
Temperature Range	+4 to 60 °C	
Temperature Deviation Over Time*	at $20 ^{\circ}\text{C} \le \pm 0.3$ at $37 ^{\circ}\text{C} \le \pm 0.3$	
Spatial Temperature Deviation*	at 20 °C: ≤ ±0.8 at 37 °C: ≤ ±0.7	at 20 °C: ≤ ±0.8 at 37 °C: ≤ ±1.0
Internal Dimensions (W x D x H)	550 x 450 x 615 mm / 21.6 x 17.7 x 24.2 in.	550 x 450 x 1015 mm / 21.6 x 17.7 x 40.0 in.
External Dimensions (W x D x H)**	652 x 605 x 1050** mm / 25.7 x 23.8 x 41.3** in.	652 x 605 x 1450** mm / 25.7 x 23.8 x 57.1** in.
Footprint	0.39 m ² / 4.2 ft ²	
Shelves Supplied / Max	2 / 4	3/5
Shelf Size (W x D)	545 x 445 mm / 21.46 x 17.52 in.	
Loading Capacity Per Shelf	15 kg / 33 lbs	
Weight Of Unit	95 kg / 209 lbs	115 kg / 254 lbs
Shipping Weight	115 kg / 254 lbs	135 kg / 298 lbs
Electrical	230 V / 50 Hz, China Plug	
Rated Power / Max. Current	900 W / 3.9 A	

^{*} Based on measurement according to DIN12880, at ambient temperature of 22 °C ± 3 °C ** Height includes the height of feet / casters

Accessories

Shelving		
Order Number	RI Shelf	
Description	Perforated stainless steel shelf (white, powder coated) with 4 shelf holders	
Dimensions	545 x 445 mm / 21.46 x 17.52 in.	
Shipping Weight	2 kg / 4 lbs	
For Use With	RI-150CN, RI-250CN	





Shelf holders included

thermofisher.com

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Australia +61 39757 4300 Austria +43 1 801 40 0 Belgium +32 53 73 42 41 China +800 810 5118 or +400 650 5118 France +33 2 2803 2180

Germany national toll free 0800 1 536 376 Germany international +49 6184 90 6000

India toll free 1800 22 8374 India +91 22 6716 2200 Italy +39 02 95059 552 Japan +81 3 5826 1616 **Netherlands** +31 76 579 55 55 **New Zealand** +64 9 980 6700 Nordic/Baltic/CIS countries +358 10 329 2200

Russia +7 812 703 42 15 **Spain/Portugal** +34 93 223 09 18 **Switzerland** +41 44 454 12 12 **UK/Ireland** +44 870 609 9203 USA/Canada +1 866 984 3766

Other Asian countries +852 2885 4613 **Countries not listed** +49 6184 90 6000



A Thermo Fisher Scientific Brand